

Where Can I Find The Installation Log-File

All the logs you see during the installation or update in the installer window are also written into a file. These files can be found in the GECAMed folder.

Installation Log-File

The installation log-file has the form *"Install_V«GECAMed version»_«YearMonthDayHourMinuteSeconds»_«a unique number».log"* (e.g. *"Install_V1.2.07_20130213100050_1234567890123456789.log"*) and lies directly in the GECAMed install folder.

Updater Log-File

There are two update log-files. The first one has the same form, as the install log file. The second is for the SQL log only. It has the form *"sql_update_«Year»_«Month»_«Day»_«Hour»_«Minute».log"*. Both files lie in the folder *"«GECAMed install folder»/updater/«GECAMed version»/"*

From:

<http://gecamed.lu/> - **GECAMed - Gestion de Cabinets Médicaux**

Permanent link:

http://gecamed.lu/faq/installation/installation_logfile

Last update: **2019/12/09 11:19**

